

PCP[®] bath safety Shower chairs



- Seats and back rests made of blow molded high impact polyethylene, completely closed
- Rustproof anodized aluminum framework, 1" (2.5 cm) o.d., with sturdy cross-x design to prevent wobbling
- Seat and backrest 12" x 20" (31 x 51 cm)
- Seat height adjustment: 16"-20" (41-53 cm)
- Backrest height from seat on models 7003, 7004-B, 7006-B: 16" (41cm)
- Shipped disassembled, 1 per carton

No. 7002 – Shower chair without back, retail package

No. 7003 – Shower chair with back, retail package

No. 7004 – Shower chair without back, heavy suction tips

No. 7004-B – Shower chair with back, heavy suction tips

No. 7006 – Shower chair without back, 15" (38 cm) seat height

No. 7006-B – Shower chair with back, 15" (38 cm) seat height

Models 7004 and 7004-B feature heavy duty 3" (7.5 cm) diameter rubber suction tips for excellent gripping action.



Transfer Bench

No. 7391

- Molded plastic seat and back panels
- Seat area 26" x 16" (66 x 41 cm) with drain holes
- Both contoured back rest (17" x 7" – 43 x 18 cm) and arm rest are reversible for right or left side
- Seat height adjusts 19"-23" (48-58 cm)
- Lightweight 1" (2.54 cm) aluminum frame
- Reinforced suction tips
- 1 per carton, 12 lbs.



Bathtub Stool

No. 7001

- Sturdy molded polystyrene seat, 14" (36 cm) diameter
- Lightweight aluminum framework
- Seat height adjusts 15"-20" (38-51 cm)
- Reinforced rubber tips
- 2 per carton, 7 lbs.



Bathtub Safety Board

No. 7397

- Makes sponge bathing easier and safer
- Injection molded board – 27" x 10" (69 x 25 cm) – with non-slip finish
- Four adjustment knobs enable bath board to fit different tub wall gaps, from 15½" (39 cm) minimum to 25" (63 cm) maximum
- Includes molded soap tray
- Drain holes for extra safety
- Side hand grips for extra support
- Non-marring tub pads



Hygienic Shower Chairs

- Seats and back rests made of blow molded high impact polyethylene, completely closed
- Large 16" x 15" (41 x 38 cm) seat with generous cutaway to allow for easy cleaning of the perineal area
- Seat height adjusts 15"-19" (38-48 cm)
- Backrest height from seat on model 7005-B: 13" (33 cm)
- Shipped disassembled, 1 per carton

No. 7005 – Shower chair without back

No. 7005-B – Shower chair with back



Cast Brackets

- Thick wall chrome steel tubing, 1" (2.5 cm) o.d., with knurled finish for secure gripping
- Brackets of heavy zinc alloy casting
- #14 x 1¾" (4 cm) long chrome plated mounting screws included with each unit


- No. 8112-K – 12" (30 cm) long
- No. 8116-K – 16" (41 cm) long
- No. 8118-K – 18" (46 cm) long
- No. 8124-K – 24" (61 cm) long
- No. 8130-K – 30" (76 cm) long
- No. 8132-K – 32" (81 cm) long
- No. 8136-K – 36" (91 cm) long
- No. 8140-LK 16" x 24" (41 x 61 cm) left angle
- No. 8140-RK – 16" x 24" (41 x 61 cm) right angle
- No. 8142-K – 42" (107 cm) long
- No. 8145-K – 12" x 12" (30 x 30 cm) 45° angle
- No. 8148-LK – 16" x 32" (41 x 81 cm) left angle
- No. 8148-RK – 16" x 32" (41 x 81 cm) right angle



Formed Brackets

- Thick wall chrome steel tubing 1" (2.5 cm) o.d., with knurled finish for secure gripping
- Strong formed steel brackets
- Mounting hardware included

- No. 9112-K – 12" (30 cm) long
- No. 9116-K – 16" (41 cm) long
- No. 9118-K – 18" (46 cm) long
- No. 9124-K – 24" (61 cm) long
- No. 9132-K – 32" (81 cm) long



Molded Plastic

- Ergonomic triangular design
- Warm to the hand
- White plastic wrapped over heavy duty aluminum tubing
- Non-slippery coating provides secure safety grip
- Heavy duty molded brackets

- No. 9212 – 12" (30 cm) long
- No. 9216 – 16" (41 cm) long
- No. 9218 – 18" (46 cm) long
- No. 9224 – 24" (61 cm) long
- No. 9232 – 32" (81 cm) long



Molded Tub Bar

No. 7250

- Extends 12" (31 cm) above tub wall
- Ergonomic design with textured surface provides for comfortable, secure gripping
- Strong molded construction
- One piece steel platform with easy cam action locking mechanism
- Wipes clean easily
- Non-marring rubber pads

Chrome Tub Bars

- Easy to install, fits all tubs with 3 1/2"-6 1/2" (9-16 cm) wall thickness
- Thick wall chrome steel tubing, 1" (2.5 cm) o.d., knurled finish for secure gripping
- Hardened die cast brackets with built-in rocker to provide self-alignment for uniform clamping surface

No. 7110-K – 6" x 11" (15 x 28 cm)
 No. 7120-K – 6" x 18" (15 x 46 cm)
 No. 7140-K – 6" & 18" x 18" (15 & 46 x 46 cm)



Coated Tub Bar

No. 7255

- Extends 15" (38 cm) above tub wall
- Strong welded construction
- One piece steel platform with stainless steel locking mechanism
- Attractive, easy to clean vinyl coating
- Non-marring rubber pads

Coated Tub Bars

- Coated thick wall steel tubing with a special finish to provide a secure grip when the hands are wet
- Strong cast metal clamps with rust resistant coating, well padded to eliminate marring of tub finish

No. 7210 – 6" x 11" (15 x 28 cm)
 No. 7220 – 6" x 18" (15 x 46 cm)
 No. 7240 – 6" & 18" x 18" (15 & 46 x 46 cm)
 No. 7260 – space saver design



3 in 1 Aluminum Commode

No. 7025

- Anodized 1" (2.5 cm) frame, 1¹/₈" (3 cm) legs
- Metal reinforced rubber tips
- Includes molded plastic seat ring, removable splash guard and no. 5114 removable 6 qt. pail with lid
- Seat height adjusts 19"-23" (48-58 cm)
- Overall width 20" (51 cm)
- Overall depth 19" (48 cm)
- Width between arms 18" (46 cm).
- 1 per carton, 12 lbs.



Safety Frame/Raised Seat

No. 7007

- Cost effective two-in-one product
- Anodized 1" (2.5 cm) frame, 1¹/₈" (3 cm) legs
- Molded plastic seat fixed to frame
- Vinyl hand grips, silencer rings, plastic splash guard
- Seat height adjusts 18"-23" (46-58 cm)
- Overall width 21" (53 cm)
- Width between arms 18" (46 cm)

Standard Safety Frame

No. 7009

- Anodized aluminum, 1" (2.5 cm) o.d.
- Full length plastic waterfall arm rests
- Locking bracket attaches with basic seat assembly bolts
- Can be lifted off floor for cleaning.
- Arm height adjusts 27"-31" (69-79 cm); width between arms adjusts 17"-22" (43-60 cm)
- Reinforced rubber tips

Deluxe Safety Frame

No. 7010

- Stable four point floor contact
- Anodized aluminum, 1" (2.5 cm) o.d.
- Adjusts for seat heights 12"-16" (30-40 cm)
- Full length plastic waterfall arm rests
- Arm height from seat 25" (63 cm); width between arms 21" (53 cm)
- Reinforced rubber tips





Raised Toilet Seat

No. 7013

- Easy to clean unbreakable polyethylene
- Wide seating area, raised in back for comfort
- Adds 5" (7.5 cm) to height of toilet seat
- Discreet carrying bag for use when traveling
- Locking mechanism not included on this model



Raised Toilet Seat With Arms

No. 7017

- Easy hand tightening mechanism safely secures seat to toilet bowl
- Lightweight and unbreakable; specially molded for sitting comfort
- Raises seat height 5" (12.5 cm)
- Lightweight, padded aluminum arm rests for extra security
- Width between arms 18" (46cm)
- Easily wipes clean
- 1 per carton, 8 lbs.



Locking Raised Toilet Seat

No. 7015

- As 7017, but does not include arm rests
- 1 per carton, 5 lbs.



Raised Toilet Seat

No. 7014

- Adjustable 3"-6" (8-15 cm)
- Easy to clean molded polyethylene ring
- Continuous polyethylene splash guard
- Aluminum legs with vinyl tips

Safety Grip Guard Seat

No. 7751

- For use in smaller bathrooms where space is at a premium
- Padded arm rests, safety clamp legs and universal legs on select models
- Other features as 7014 at left



7751



Hand-Held Shower, Standard
No. 7031

- Attaches easily to standard 1/2" shower arm
- Simple on/off valve on water spray controls flow of water
- Nylon reinforced hose allows water spray to be used while sitting on shower seat or bench
- Wall mounting bracket can be positioned anywhere on wall
- Compatible with no. 7034 diverter valve

Hand-Held Shower, Deluxe
No. 7033

- Water spray features easy on/off valve with adjustable flow control
- Reinforced vinyl hose allows spray to reach all parts of user's body while bathing
- Attaches easily to standard 1/2" shower arm
- Wall mounting bracket included
- Compatible with no. 7034 diverter valve

Shower Diverter Valve
No. 7034

- Attaches quickly and easily to shower arm; a turn of the valve allows use of either regular shower head or hand held shower
- Compatible with either 7031 or 7033



Commode Cushion
No. 7018

- Comfortable foam padding, 2" thick, covered with naugahyde – wipes clean easily
- Provides extra comfort and helps minimize pressure sores
- Velcro fastening straps allow for easy attachment to and removal from toilet seat
- Fits all standard toilet seats



Bath Mats

No. 7036 – Ivory
No. 7037 – Mauve
No. 7038 – Blue

- Extra large suction cups keeps mat from slipping on tub floor
- Extra large size – 16" x 31" – covers most of tub floor
- Rubber non-slip mat stays soft and pliant
- Mildew resistant
- Machine washable