



**Gray Coated 3-in-1 Commode**

- 7/8" (2 cm) bolted steel frame, 1" (2.5 cm) legs
- Easily removable angled back rest
- Molded plastic arm rests for comfort
- Metal reinforced rubber tips
- Snap-on seat assembly and removable half-pail with lid
- Seat height adjustment: 16"-22" (41-56 cm)
- Overall width 26" (66 cm)
- Overall depth 20" (51 cm)
- Width between arms 18" (46 cm)
- 1 per carton, 18 lbs.

No. 5526-ESU – Shipped assembled, 1 per carton

No. 5526-RBE – Shipped disassembled, 2 per carton



**Deluxe Aluminum Commode**

No. 5027

- Angled back, plastic waterfall arm rests for greater comfort
- Anodized 1" (2.5 cm) frame, 1 1/8" (3 cm) legs
- Metal reinforced rubber tips
- Includes no. 5116 seat assembly and no. 5114 removable 6 qt. pail with lid
- Seat height adjustment: 19"-23" (48-58 cm)
- Overall width 22" (56 cm)
- Overall depth 20" (51 cm)
- Width between arms 17" (43 cm)
- 1 per carton, 17 lbs.



**Removable Back Commode**

No. 5026

- Easily removable angled back rest
- Molded plastic arm rests for comfort
- Anodized 1" (2.5 cm) frame, 1 1/8" (2.5 cm) legs
- Metal reinforced rubber tips
- Includes no. 5116 seat assembly and no. 5114 removable 6 qt. pail with lid
- Seat height adjustment: 19"-23" (48-58 cm)
- Overall width 21" (53 cm)
- Overall depth 20" (51 cm)
- Width between arms 17" (43 cm)
- 1 per carton, 14 lbs.



**Aluminum Commode**

No. 5018

- Anodized 1" (2.5 cm) frame, 1 1/8" (3 cm) legs
- Metal reinforced rubber tips
- Includes no. 5116 seat assembly and no. 5114 removable 6 qt. pail with lid
- Seat height adjustment: 19"-23" (48-58 cm)
- Overall width 22" (56 cm)
- Overall depth 20" (51 cm)
- Width between arms 18" (46 cm)
- 1 per carton, 14 lbs.



5028

### Drop Arm Commode

No. 5028

- Arms drop independently below seat level for safe and easy lateral transfers
- Convenient arm release locks automatically in raised position
- Molded plastic arm rests for comfort
- $\frac{7}{8}$ " (2 cm) welded frame, 1" (2.5 cm) legs
- Metal reinforced rubber tips
- Includes no. 5116 seat assembly and no. 5114 removable 6 qt. pail with lid
- Seat height adjustment: 19"-23" (48-58 cm)
- Overall width 21" (53 cm)
- Overall depth 19" (48 cm)
- Width between arms 17" (43 cm)
- 1 per carton, 18 lbs.



5029

### Extra-Wide Commode

No. 5029

- $\frac{7}{8}$ " (2 cm) welded chrome frame, 1" (2.5 cm) legs
- Molded plastic arm rests for comfort
- Metal reinforced rubber tips
- Includes no. 5116 seat assembly and no. 5114 removable 6 qt. pail with lid
- Seat height adjustment: 19"-22" (48-56 cm)
- Overall width 28" (71 cm)
- Overall depth 21" (53 cm)
- Width between arms 24" (61 cm)
- 1 per carton, 21 lbs.



5004

### Mobile Shower Chair

No. 5004

- Allows for virtually total access to user's body during bathing
- Flexible nylon backrest
- Anodized 1" (2.5 cm) frame,  $1\frac{1}{8}$ " (2.5 cm) legs
- Rust resistant 4" (10 cm) ball bearing casters with rear wheel locks
- Includes no. 5116 seat assembly and no. 5114 removable 6 qt. pail with lid
- Once pail is removed, chair moves easily over most toilets
- Seat height 21" (53 cm)
- Overall width 21" (53 cm)
- Overall depth 20" (51 cm)
- Width between arms 18" (46 cm)
- 1 per carton, 21 lbs.



### Pail with Lid and Handle

No. 5115

- 10-quart capacity



### Pail with Lid

No. 5114

- 6-quart capacity



### Seat Assembly

No. 5116

- Snaps on and off for easy cleaning



Quick and easy set up: no bolts on the cross bars or brackets.  
All attachments made by spring loaded clips and buttons.



Pivot drop feature allows one end to be lowered while the other remains in the up position.



### Telescopic Bed Rails

No. 6081

- Designed for quick and easy installation on any hospital bed frame with width from 24" to 44" (61-112 cm)
- Set telescopes from full length to three quarter length, 48"-74" (122-188 cm)
- Chrome steel 1" (2.5 cm) tubing
- Clamps feature spring loaded clips and buttons
- Height above spring, up: 16" (41 cm)
- Height above spring, down: 1" (2.5 cm)
- Net weight of set: 28 lbs.
- Shipped 1 set per carton, 33 lbs.

### Half-Length Rails

No. 6083

- Fits all standard hospital beds – no tools required for installation
- Anodized aluminum 1" (2.5 cm) tubing
- Overall length, 37" (94 cm)
- Height above spring, up: 16" (41 cm)
- Net weight of set: 13 lbs.
- Shipped 1 set per carton, 14 lbs.



### Contour Style Mattress Cover

No. 6211

- Fits hospital size bed

### Zipper Style Mattress Cover

No. 6213

- Covers entire hospital bed size mattress



**Sitz Bath**

No. 5040

- Vinyl solution bag can be hung over commode or set on top of tank
- Molded plastic bowl with wide contour edge and tubing anchor at bottom
- Compatible with most commodes

**Male Urinal**

No. 6091

- Molded plastic, one quart capacity
- Hinged cover, open handle

**Fracture Pan**

No. 6095

- Molded plastic with loop handle
- For use with fracture or immobile patients, also as female urinal

**Bed Pan**

No. 6094

- Molded plastic
- Wide flange for easy gripping

**Foot Stool**

No. 5008

- Sturdy and tip proof
- Chrome steel legs, metal reinforced rubber tips
- Rugged non-skid 10" x 14" (25 x 36 cm) top
- Height 9" (23 cm)

**Blanket Support**

No. 6210

- Keeps blanket off feet to protect sensitive skin and prevent chapping
- Wire steel frame, adjusts up to 19" (48 cm), drops to top of mattress when not in use
- Base inserts between spring and mattress

**Shampoo Tray**

No. 6096

- Inflatable ring allows for easy shampooing of bed confined patient
- Water carriage tube minimizes spills


**Waterproof Sheeting**

- Reversible, washable, moisture proof

No. 6258 – 36" x 54" (91 x 137 cm)

No. 6259 – 36" x 66" (91 x 168 cm)






### Bed Wedges

- Provides gradually sloping support for head, foot or leg elevation
- Firm density foam with white dacron/cotton cover, removable for washing
- Completely non-allergenic
- Distinctive consumer packaging

No. 6157 – 7<sup>1</sup>/<sub>2</sub>" x 26" x 24"  
(19 x 66 x 61 cm)

No. 6158 – 10<sup>1</sup>/<sub>2</sub>" x 26" x 24"  
(27 x 66 x 61 cm)

No. 6159 – 12<sup>1</sup>/<sub>2</sub>" x 26" x 24"  
(32 x 66 x 61 cm)



### Deluxe Pressure Relief Bed Pads

- Stain resistant natural color polyester/acrylic, will not pack down
- Stays soft and thick through repeated washing

No. 6261 – 24" x 30" (61 x 76 cm)

No. 6262 – 30" x 40" (76 x 102 cm)

No. 6263 – 30" x 60" (76 x 152 cm)

No. 6264 – Hospital Bed Size

No. 6266 – Wheelchair Pad

### Standard Bed Pads

- 100% white polyester
- Totally machine washable

No. 6203 – 24" x 30" (61 x 76 cm)

No. 6204 – 30" x 40" (76 x 102 cm)

No. 6205 – 30" x 60" (76 x 152 cm)

No. 6206 – Heel Protectors (1 pr)

No. 6207 – Elbow Protectors (1 pr)

No. 6208 – Wheelchair Pad

### Egg Crate Bed Pads

- Convoluted foam pad, 25" x 74" (89 x 190 cm)
- Helps minimize pressure sores by equalizing weight distribution

No. 6222 – 2" (5 cm) thick

No. 6223 – 3" (7.5 cm) thick