



Three Wheeled Rollator

No. 5303-BL – Blue Finish

No. 5303-CH – Chrome Finish

- Large 8" pneumatic wheels
- Anatomical hand grips, adjustable in height from 33" to 38" with quick lock handles
- Loop-type cable hand brakes operate on rear wheels
- Includes wire basket, 9" x 9" x 3", easily detachable from frame
- Easy folding mechanism for transport and storage
- 14 lbs. total



Four Wheeled Rollator with Curved Back Rest

No. 5313-TI – Titanium Finish

- Curved backrest for more comfort when sitting
- Comfortably padded seat
- Basket can be positioned either under or in front of seat
- Adjustable handle height for custom fit
- Easy to operate ergonomic hand brakes lock rear wheels for security
- Truly effortless folding mechanism
- Large 8" pneumatic wheels for smoother rolling operation
- 17" between hand grips
- Product weight 19 lbs



Four Wheeled Rollator with Back Rest

No. 5314-BG – Burgundy Finish

No. 5314-GR – Teal Green Finish

- Adjustable height handles with foam hand grips: adjust 34"-38"
- Loop type cable hand brakes operate on rear wheels
- Padded seat lifts up to accommodate removable 6" x 6" x 13" wire basket with handle
- Smooth rolling 5" wheels
- Padded removable backrest for extra comfort
- Easy folding mechanism for transport and storage
- Lightweight – 12 lbs. total



5066-E

- Folds easily for storage and transport
- Palm activated top release folding mechanism
- 1" (2.5 cm) frame, 1 1/8" (3 cm) legs
- Steel front cross brace with internal bushing for smooth operation
- Flush side cross braces for added stability
- Metal reinforced rubber tips.
- Adult and junior (add suffix "-J") sizes



5060-05

- No. 5066-E – Coated front cross brace
- No. 5066-TR – Chrome front cross brace
- No. 5566-B – Bronze frame, texture coated front cross brace
- No. 5060-05 – Coated front cross brace, 5" swivel wheel attachments with rear auto glide legs
- No. 5062-TR – Coated front cross brace, 3" single fixed wheel attachments in front



5069

Extra Wide Folding Walker

No. 5069

- Allows bariatric user to manipulate walker safely and comfortably
- Width between hand grips 19 1/2" (50 cm), 2 1/2" (6 cm) wider than standard
- Bottom safety release folding mechanism
- All other design features and materials as 5066-TR
- Adult size only



5566-B



5062-TR



- Folds easily for storage and transport
- Palm activated push button folding mechanism
- 1" (2.5 cm) frame, 1 1/8" (3 cm) legs
- Steel front cross brace with internal bushing for smooth operation
- Double side flush cross braces for optimum stability
- Metal reinforced rubber tips
- Adult and junior (add suffix "-J") sizes – model 5050 only



No. 5050-E – Coated front cross brace
 No. 5052-E – Coated front cross brace, 5" single fixed wheel attachments in front



Walker Ski Glides

No. 5106

- Allows user to maneuver basic obstacles like uneven sidewalks and other challenging surfaces.
- High density plastic with a nylon surface.
- Sized to fit 1", 7/8", 3/4", with 5/8" o.d.

Lightweight Walker

No. 5055-E

- Considerably lighter in weight, yet as stable as standard walkers
- 7/8" (2.2 cm) frame, 3/4" (2 cm) legs
- Coated front cross brace
- Adult size only
- 1 per carton, 6 lbs.
- All other features as model 5050




L to R: 5059, 5059-E, 5109-E

Wheel Attachments

No. 5059 – 3" (7.5 cm) fixed wheel attachments, 1" (2.5 cm) wide

No. 5059-E – 3" (7.5 cm) fixed wheel attachments, 1/2" (1.3 cm) wide

No. 5109-E – 5" (12 cm) fixed wheel attachments, fits models 5050-E, 5066, 5069



L to R: 5107, 5103, 5105

Wheel Attachments

No. 5107 – 5" (12 cm) swivel wheel attachments

No. 5103 – 3" (7.5 cm) swivel wheel attachments

No. 5105 – Auto glide attachment, nylon tip spring loaded guide riveted onto leg piece



Rigid Adjustable Walker
No. 5067

- Inverted U-design front legs for strength and stability
- Angled legs provide for wide center of gravity, added stability and tipping resistance
- Internal reinforcement at all stress points
- Anodized 1" (2.5 cm) frame, 1 1/8" (3 cm) legs
- Cushioned hand grips
- Adult, junior (5067-J) and child (5067-C) sizes
- Shipped 2 per carton



Folding Hemi Walker
No. 5063

- For users with limited or no dexterity in one arm or hand
- 7/8" (2.2 cm) frame
- Folding cross brace with firm locking action
- Angled legs for full tip contact with floor
- Black hypalon hand grips
- Adult size only



Walker Platform Attachment
No. 5049

- Indicated when weight cannot be borne by the wrist
- Attaches easily to walker by means of two horizontal side tightening knobs; no tools required
- Easy height adjustment; reversible for right or left hand use
- Can be mounted on inside or outside of frame
- Padded vinyl cuff with hook and loop fastening
- Handle can be rotated to desired angle
- Lightweight aluminum and steel construction



Walker Pouch
No. 5164

- Blue denim fabric with hook and loop fastening
- Large 10" x 15" (25 x 38 cm) pocket with smaller dual inner pockets.
- Machine washable



Walker Basket
No. 5166

- Fits model nos. 5050, 5052 only
- White vinyl coated steel wire
- Molded plastic tray at bottom



L to R: 5073, 5092, 5090

Wood Crutches

- No. 5073 – Adult, for heights 5'4"-6'5"
- No. 5074 – Medium, for heights 4'1"-5'7"
- No. 5075 – Junior, for heights 4'5"-5'
- No. 5076 – Child, for heights 3'8"-4'4"
- Includes tips, grips, and axillary pads shrink wrapped with user instruction card

Aluminum Crutches

- No. 5091 – Tall, for heights 5'10"-6'6"
- No. 5092 – Adult, for heights 5'2"-5'10"
- No. 5093 – Junior, for heights 4'6"-5'2"
- No. 5094 – Child, for heights 4'2"-4'8"
- Includes tips, grips, and axillary pads shrink wrapped with user instruction card
- Easy push button adjustment with user height scale

Forearm Crutches

- No. 5090 – Adult, for heights 5'2"-6'2"
- No. 5090-J – Junior, for heights 4'2"-5'2"
- 1" (2 cm) o.d. aluminum tubing.
- Coated cuffs with 180° pivot
- Molded black vinyl hand grips on welded and coated steel handles
- Adjustment at legs and forearm cuffs



Pyramid Base Quad Canes

- No. 5083 – Three point (8" dia.) base
- No. 5084 – Large (9" x 12") quad base
- No. 5087 – Small (6" x 8") quad base
- Chrome steel 3/4" (2 cm) tubing
- Bolt and wing nut adjustment
- Offset handle, gray vinyl handgrip
- Height adjustment: 30"-36" (76-91 cm).
- Shipped 2 per carton



Low Profile Quad Canes

- 7/8" (2.2 cm) aluminum shaft
- Steel base: small (6" x 8", 15 x 20 cm) or large (8" x 12", 20 x 30 cm)
- Offset handle, molded gray vinyl or black hypalon foam handgrip
- Safety lock silencer ring
- Shipped 2 per carton

- No. 5095-E – Small base, silver frost, gray vinyl grip
- No. 5096-E – Large base, silver frost, gray vinyl grip
- No. 5595-BL – Small base, black
- No. 5596-BL – Large base, black
- No. 5595-BP – Small base, blue pewter
- No. 5596-BP – Large base, blue pewter
- No. 5595-BR – Small base, bronze
- No. 5596-BR – Large base, bronze

**Proudly introducing a completely retooled line of mobility aids
 – most MADE-IN-THE-U.S.A. – priced to compete with the imports!**



A. Offset handle with lanyard:

- No. 5182 – paisley
- No. 5183 – compass rose
- No. 5184 – butterfly

B. Offset handle:

- No. 212400 – silver frost
- No. 212200 – blue pewter
- No. 212100 – black

C. Medium crook handle:

- No. 211400 – silver frost
- No. 211800 – bronze
- No. 211600 – blue pewter

D. Large crook handle:


- No. 210600 – silver frost
- No. 210300 – bronze
- No. 212300 – blue pewter

E. Medium orthopedic handle:

- No. 212000 – silver frost
- No. 212800 – bronze
- No. 210500 – blue pewter

F. Large orthopedic handle:

- No. 211500 – silver frost
- No. 210800 – bronze
- No. 210200 – blue pewter
- No. 210100 – black



Available colors

-  • silver frost
-  • blue pewter
-  • black
-  • bronze
-  • paisley
-  • compass rose
-  • butterfly



Adjustable folding canes

Large orthopedic handle:

- No. 210400 – bronze
- No. 210402 – blue pewter
- No. 210404 – black

Medium orthopedic handle:

- No. 210403 – bronze
- No. 210411 – blue pewter
- No. 210412 – black



Highly visible and easily accessible Mini-Departments (Nos. 4071, 4072)

attach to either slatwall or perfbboard, and feature an ideal selection of canes. The display unit keeps merchandise organized for easy purchase, and our unique labeling system allows quick, easy and accurate restock.

- Ergonomically designed hand grips
- Gender specific sized handles (medium and large)
- Available in attractive, non-chip colors and patterns that combine form with function
- Safety locking bolt on quad canes
- Convenient contact closure fastening strap on folding styles eliminates the need for carrying case
- Rubber friction tip on orthopedic handles
- First-quality tips




A. Medium T-handle:

No. 5122 – black handle, walnut
 No. 5146 – brown handle, walnut

B. Large T-handle:

No. 5123 – black handle, walnut
 No. 5148 – brown handle, walnut




C. Medium derby handle:

No. 5131-02 – light brown
 No. 5131-12 – ramen wood
 No. 5131-20 – beechwood
 No. 5131-32 – cherrywood

D. Large derby handle:

No. 5131-01 – light brown
 No. 5131-11 – ramen wood
 No. 5131-21 – beechwood, black



G. Medium crook handle:

No. 5120 – walnut (7/8")
 No. 5121 – mahogany (7/8")
 No. 5135 – walnut
 No. 5140 – rosewood

H. Large crook handle:

No. 5126 – walnut
 No. 5127 – mahogany
 No. 5133-41 – stained chestnut
 No. 5160 – carved veneer



PCP's Mobility Central

Merchandiser No. 4075 requires only 2.5 square feet of retail space, holds three different categories of product, including a wide variety of canes as well as the most commonly requested accessory items. Tearaway usage instructions both educates and sells.




E. Medium Fritz handle:

No. 5132-02 – light brown
 No. 5132-22 – beechwood, black

F. Large Fritz handle:

No. 5132-01 – light brown
 No. 5132-21 – beechwood, black



6107 Ice Cane Attachment

- For safer walking in ice and snow.
- Strong steel five prong tip.

Available colors

- beechwood, black
- cherrywood
- chestnut
- light brown
- mahogany
- ramen wood
- rosewood
- walnut



Tips and Accessories

- No. 5070 crutch pads
- No. 5071 crutch grips, open style
- No. 5072 crutch grips, closed style
- No. 6099 crutch tips, small
- No. 6100 crutch tips, medium
- No. 6101 crutch tips, large
- No. 6102 Walker tips, fits 1 1/8"
- No. 6103 Cane tips, fits 1"
- No. 6104 Cane tips, fits 7/8"
- No. 6105 Cane tips, fits 3/4"
- No. 6106 Quad cane tips, fits 5/8"