

CERVICAL COLLARS



MODEL 9977

Soft foam collar for mild cervical lesions, or wherever a reminder against making quick and painful neck movements is required. Washable polyurethane foam covered with ribbed tubular knit stockinette; hook-and-loop fastening. Each size available in three depths: narrow (2.5"), average (3"), wide (3.5").

Sizes – X-Small: 11.5"-13.5" neck
Small: 13.5"-15.5" neck
Medium: 15.5"-17.5" neck
Large: 17.5"-19.5" neck
X-Large: 19.5"-21.5" neck



MODEL 9976

Firm, washable polyurethane foam covered with beige two-way stretch hosiery stockinette. Obscures soiling, enhances durability for effective longer-term wear. Hook-and-loop fastening. Same sizes as model 9977 above.



MODEL 975

For flexion and hyperextension (whiplash), injuries, torticollis and osteoarthritis. Perforated polyethylene with vinyl covered foam padded edges; hook-and-loop fastening and adjustment panels. Simple and quick anterior height adjustment – from 3" to 5". Adapts for flexion, extension and torticollis. Same sizes as model 9977 at left.



MODEL 976

Thomas type collar; same general design and components as 975 above, but with positive screw-lock height adjustment slots – adjusts 3" to 5". Same sizes as model 9977 at left.



MODEL 9979

For moderate cervical lesions, including whiplash injuries, osteoarthritis, and cases of mild to moderate soft tissue damage. Perforated polyethylene with vinyl covered foam padded edges; hook-and-loop fastening and adjustment panels. Contoured chin rest with removable pad enhances wearer comfort, provides effective flexion control. Anterior height adjusts from 4.5" to 6" – sternal notch to lower mandible. Same sizes as model 9977 at left.



MODEL 977

Same design and application as 975, but made from cool, breathable, semi-rigid mesh material – ideal for wear in warmer climates. Hook and loop anterior height adjustment – adapts for flexion, extension and torticollis.

