

MEN'S LUMBOSACRAL SUPPORTS, LACED



MODEL 464-GHS

Widely used, basic support with broad range of applications for the slighter-than-average figure.



MODEL 464-GE

Same design as 464-GHS, but made from gray material to help obscure soiling. Appropriate for the more active wearer.



MODEL 466-HS

Designed for the short, stocky figure – 3 side pulls to enhance cosmesis, 4 steels in back for more effective motion control.



MODEL 464-VHS

Same design and materials as 464-GHS, but with contact closure tabs in front for the wearer with limited dexterity or hand strength.



MODEL 464-SHS

Similar in design to 464-GHS, but shorter both in front and back. Appropriate for the adolescent or someone who will be primarily sitting while wearing the corset.



MODEL 829-HS

Similar in design to 464-GHS, but deeper both in front and back for the slightly taller, wiry figure type. Also recommended where inventory control is a top priority.

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	ADJUSTMENT	FRONT CLOSURE	ELASTIC RELEASES	STEELS	EVEN HIP SIZES
464-GHS	Thin to average	Short	FRONT 7" BACK 12"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 10"	28-50
464-VHS	Thin to average	Short	FRONT 7" BACK 12"	White canvas	Two pull side lace	Hook-and-loop contact tabs	At groin	1 pr. 10"	28-50
464-GE	Thin to average	Short	FRONT 7" BACK 12"	Gray canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 10"	28-50
464-SHS	Thin to average	Short	FRONT 7" BACK 10"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 8"	28-44
466-HS	Stocky	Short	FRONT 8" BACK 12"	White canvas	Three pull side lace	Snap button/ hook-and-eye	At top and groin	2 pr. 10"	30-48
829-HS	Thin to average	Average	FRONT 9" BACK 15"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 13"	28-48

MEN'S LUMBOSACRAL SUPPORTS, LACED



MODEL 464-GHS

Widely used, basic support with broad range of applications for the slighter-than-average figure.



MODEL 464-GE

Same design as 464-GHS, but made from gray material to help obscure soiling. Appropriate for the more active wearer.



MODEL 466-HS

Designed for the short, stocky figure – 3 side pulls to enhance cosmesis, 4 steels in back for more effective motion control.



MODEL 464-VHS

Same design and materials as 464-GHS, but with contact closure tabs in front for the wearer with limited dexterity or hand strength.



MODEL 464-SHS

Similar in design to 464-GHS, but shorter both in front and back. Appropriate for the adolescent or someone who will be primarily sitting while wearing the corset.



MODEL 829-HS

Similar in design to 464-GHS, but deeper both in front and back for the slightly taller, wiry figure type. Also recommended where inventory control is a top priority.

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	ADJUSTMENT	FRONT CLOSURE	ELASTIC RELEASES	STEELS	EVEN HIP SIZES
464-GHS	Thin to average	Short	FRONT 7" BACK 12"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 10"	28-50
464-VHS	Thin to average	Short	FRONT 7" BACK 12"	White canvas	Two pull side lace	Hook-and-loop contact tabs	At groin	1 pr. 10"	28-50
464-GE	Thin to average	Short	FRONT 7" BACK 12"	Gray canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 10"	28-50
464-SHS	Thin to average	Short	FRONT 7" BACK 10"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 8"	28-44
466-HS	Stocky	Short	FRONT 8" BACK 12"	White canvas	Three pull side lace	Snap button/ hook-and-eye	At top and groin	2 pr. 10"	30-48
829-HS	Thin to average	Average	FRONT 9" BACK 15"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 13"	28-48

MEN'S LUMBOSACRAL SUPPORTS, LACED



MODEL 454-HS

Elastic side panels provide for superior cosmesis and comfort, operate independently of side pull straps. Zipper closure for easy application and removal. Casings for additional steels if desired.



MODEL 441-HS

Same design and materials as 454-HS, but with contact closure tabs in front for the wearer with limited dexterity or hand strength.



MODEL 368-E SACROILIAI SUPPORT

Suitable for virtually all figure types. Includes removable sacral pad, adjustable lace back and single encircling pull strap. Three anterior pull straps with truss buckle fastening in front. Elastic releases at groin. Gray material obscures soiling. Front depth 6", back depth 9". Even hip sizes 28-46.



MODEL 453-HS

Similar in design and same materials as 454-HS, but shorter both in front and back. Casings for additional steels if desired.



MODEL 467-HS

Hoke type garment is effective for control over muscular laxity in the torso – its broad size range allows for application even in adolescents.



MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	ADJUSTMENT	FRONT CLOSURE	ELASTIC RELEASES	STEELS	EVEN HIP SIZES
454-HS	Average to stocky	Average	FRONT 10" BACK 15"	Canvas/ elastic	Three pull side lace	Zipper	Side panels all elastic	1 pr. 13"	30-48
441-HS	Average to stocky	Average	FRONT 10" BACK 15"	Canvas/ elastic	Three pull side lace	Hook-and-loop contact tabs	Side panels all elastic	1 pr. 13"	30-48
453-HS	Thin to average	Short	FRONT 9" BACK 11"	Canvas/ elastic	Three pull side lace	Zipper	Side panels all elastic	1 pr. 9"	28-44
467-HS	Average	Average	FRONT 11.5" BACK 16.5"	White canvas	Pull strap/ truss buckle	Hoke type	At top and groin	2 pr. 13"	28-44

MEN'S DORSOLUMBAR SUPPORTS, LACED



MODEL 486-HS

Basic dorsolumbar support with broad range of application for the thin to average figure type.




MODEL 488-HS

Similar in design to 486-HS, but shorter in both front and back for the slighter than average figure type.



MODEL 487-HS

Wide top back panel to accommodate four steels if necessary. Appropriate for the heavier figure. Casings for additional steels if desired.



MODEL 591 TAYLOR (TLSO) INSERT

Converts a lumbosacral corset to a dorsolumbar by inserting the aluminum uprights into the stay casings and affixing the extension to the top of the corset (all attachments included). Removable and washable padded cover, foam shoulder straps.

Lengths: 17", 19", 21"

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	ADJUSTMENT	FRONT CLOSURE	ELASTIC RELEASES	STEELS	EVEN HIP SIZES
486-HS	Thin to average	Average	FRONT 9" BACK 24"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 20"	30-48
488-HS	Thin to average	Short	FRONT 8" BACK 15"	White canvas	Two pull side lace	Snap button/ hook-and-eye	At groin	1 pr. 18"	30-48
487-HS	Average to stocky	Average	FRONT 11.5" BACK 24"	White canvas	Three pull side lace	Snap button/ hook-and-eye	At top and groin	1 pr. 20"	30-48

SHOULDER AND ABDOMINAL SUPPORTS



MODEL 902 CLAVICLE STRAP

Indicated for clavicle fracture or wherever control of shoulder motion is required; designed for continual wear. Y-shape strap arrangement: padded stockinette covered shoulder section, non-elastic web posterior straps. Triple claw buckle posterior fastening around metal D-ring – padded back to prevent skin irritation.

Sizes – X-Small: pediatric
Small: average female
Medium: average male
Large: large male



MODEL 906 POSTURE CORRECTOR

Worn as a reminder to minimize drooping of the shoulder resulting from poor posture; when the shoulders slump the snugly adjusted axillary straps cause mild discomfort, thus encouraging correct posture. Cloth back panel reinforced with flexible boning. Adjustable web shoulder straps with axillary pads. Chest encircling strap with truss buckle fastening.

Sizes – X-Small: to 26" chest
Small: 27"-32" chest
Medium: 33"-38" chest
Large: 39"-44" chest
X-Large: over 45" chest



MODEL 414 OBESITY SUPPORT

To support the very pendulous abdomen: abdominal uplift is effected by four side front pull straps and augmented by straps crossing the hips. White canvas front and back panels, power net elastic sides. Heavy duty truss buckle fastening. 12" front, 10" back.

Even waist sizes 34-54.



MODEL 959 HERNIA SUPPORT

For use after hernia repair to prevent recurrence. Open lace back with bilateral uplift pull. Perineal straps included. Pockets for removable hernia pads (pads not included). 6.5" front, 5" back; snap button fastening reinforced with hooks and eyes.

Even waist sizes 28-44.



MODEL 426 OBESITY SUPPORT

For use after abdominal surgery, or when abdominal muscles have lost their tone. Power net elastic front and sides, white canvas back. Back lace, bilateral pull straps. 7"-9" front, 6.5"-8" back; snap button fastening reinforced with hooks and eyes.

Even hip sizes 28-44.



MODEL 958 HERNIA SUPPORT (SINGLE)

For use after hernia repair to prevent recurrence. Strong 4" elastic body band with contact closure fastening. Removable foam pad in fabric panel. Perineal strap with snap button fastening. Specify for right or left side.

Sizes – Small: 30"-32" hips
Medium: 34"-36" hips
Large: 38"-40" hips
X-Large: 42"-44" hips

MODEL 958-D HERNIA SUPPORT (DOUBLE)

As 958 above, but designed for use after bilateral hernia repair.