

## WOMEN'S LUMBOSACRAL SUPPORTS, LACED



### MODEL 613-HS

Widely used, basic support with broad range of applications for the slighter-than-average figure.



### MODEL 713-HS

Same design as 613-HS, but made from cool, lightweight mesh material – ideal for wear in warmer climates.



### MODEL 612-HS

### MODEL 612-VHS (SHOWN)

Very similar in design to 613-HS, but cut to accommodate a very straight 2" hip development. Both standard (612-HS) and contact closure tabs (612-VHS) available.



### MODEL 613-VHS

Same design and materials as 613-HS, but with contact closure tabs in front for the wearer with limited dexterity or hand strength.



### MODEL 620-HS

Very similar in design to 613-HS, but cut to accommodate very straight hips, with only a 4" development.

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	ADJUSTMENT	FRONT CLOSURE	ELASTIC RELEASES	STEELS	HIP DEV.	EVEN HIP SIZES
613-HS	Thin to average	Petite	FRONT 8" BACK 12"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 10"	6"	30-50
613-VHS	Thin to average	Petite	FRONT 8" BACK 12"	White brocade	Two pull side lace	Contact closure tabs	At groin	1 pr. 10"	6"	30-50
713-HS	Thin to average	Petite	FRONT 8" BACK 12"	White mesh	Two pull side lace	Hook-and-eye	At groin	1 pr. 10"	6"	30-50
620-HS	Straight	Petite	FRONT 6" BACK 12"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 10"	4"	30-48
612-HS	Straight	Petite	FRONT 6" BACK 12"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 10"	2"	30-50
612-VHS	Straight	Petite	FRONT 6" BACK 12"	White brocade	Two pull side lace	Contact closure tabs	At groin	1 pr. 10"	2"	30-50

## WOMEN'S LUMBOSACRAL SUPPORTS, LACED



**MODEL 616-HS**

**MODEL 616-VHS (SHOWN)**

Similar in design to 613-HS, but with a slightly longer back and more of a cupped front with elastic releases at top. Both standard (616-HS) and contact closure front (616-VHS) available.



**MODEL 628-HS**

As 616-HS, but made from cool, lightweight mesh material – ideal for wear in warmer climates.



**MODEL 619-HS**

Basic two-pull support designed for the average to tall figure with comparatively straight hips



**MODEL 617-HS**

Same design as 619-HS, but made from cool, lightweight mesh material – ideal for wear in warmer climates.



**MODEL 3619-HS**

Similar to 619-HS, but with design and material upgrades: improved cupping around buttocks, closer contouring at top, elastic binding top and bottom, tricot tongues under lacing for comfort



**MODEL 618-HS**

Very similar in design to 619-HS – slightly shorter front and back depths.

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	ADJUSTMENT	FRONT CLOSURE	ELASTIC RELEASES	STEELS	HIP DEV.	EVEN HIP SIZES
616-HS	Thin to average	Petite	FRONT 8" BACK 13"	White brocade	Two pull side lace	Hook-and-eye	At top and groin	1 pr. 11"	6"	30-44
616-VHS	Thin to average	Petite	FRONT 8" BACK 13"	White brocade	Two pull side lace	Contact closure tabs	At top and groin	1 pr. 11"	6"	30-44
628-HS	Thin to average	Petite	FRONT 8" BACK 13"	White mesh	Two pull side lace	Hook-and-eye	At top and groin	1 pr. 11"	6"	30-44
619-HS	Thin to average	Average to tall	FRONT 11" BACK 15"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 12"	6"	30-44
617-HS	Thin to average	Average to tall	FRONT 11" BACK 15"	White mesh	Two pull side lace	Hook-and-eye	At groin	1 pr. 12"	6"	30-44
3619-HS	Thin to average	Average to tall	FRONT 10" BACK 15"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 11"	6"	30-44
618-HS	Thin to average	Average to tall	FRONT 10" BACK 14.5"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 11"	6"	30-44

## WOMEN'S LUMBOSACRAL SUPPORTS, LACED



### MODEL 609-HS

Widely used, basic three-pull LS with broad range of applications for various figure types. Provides superior cosmesis and adaptability for a wide range of situations.



### MODEL 649-HS

Same design as 609-HS, but made from cool, lightweight mesh material – ideal for wear in warmer climates.



### MODEL 3611-HS

Same design and materials as 3609-HS, but includes a snap-out panty piece and wider releases both top and bottom.



### MODEL 609-VHS

Same design and materials as 609-HS, but with contact closure tabs in front for the wearer with limited dexterity or hand strength.



### MODEL 3609-HS

Similar to 619-HS, but with design and material upgrades: improved cupping around buttocks, closer contouring at top, elastic binding top and bottom, tricot tongues under lacing for comfort.



### MODEL 809-HS

Similar in design to 609-HS, but with a 2" shorter front, for more comfort when sitting.

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	ADJUSTMENT	FRONT CLOSURE	ELASTIC RELEASES	STEELS	HIP DEV.	EVEN HIP SIZES
609-HS	Average	Average	FRONT 10" BACK 16"	White brocade	Three pull side lace	Hook-and-eye	At top and groin	1 pr. 13"	8"	30-56
609-VHS	Average	Average	FRONT 10" BACK 16"	White brocade	Three pull side lace	Contact closure tabs	At top and groin	1 pr. 13"	8"	30-56
649-HS	Average	Average	FRONT 10" BACK 16"	White mesh	Three pull side lace	Hook-and-eye	At top and groin	1 pr. 13"	8"	30-48
3609-HS	Average	Average	FRONT 9.5" BACK 12"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 11"	8"	30-48
3611-HS	Average	Average	FRONT 9.5" BACK 12"	White brocade	Two pull side lace	Hook-and-eye	At groin	1 pr. 11"	8"	30-48
809-HS	Average	Average	FRONT 8" BACK 16"	White brocade	Three pull side lace	Hook-and-eye	At top and groin	1 pr. 13"	8"	30-48

## WOMEN'S LUMBOSACRAL SUPPORTS, LACED



### MODEL 615-HS

Basic three-pull LS, like 609-HS, but with a cupped front to provide more abdominal uplift for the slightly stockier figure.



### MODEL 702-HS

A groin length garment which features a generous elastic gusset at the bottom front. Good for a heavier figure that does not require all the coverage of a girder length garment.



### MODEL 639-HS

Same design and materials as 630-HS, but with contact closure tabs in front for the wearer with limited dexterity or hand strength. Casings for additional steels if desired.



### MODEL 615-VHS

Same design and materials as 613-HS, but with contact closure tabs in front for the wearer with limited dexterity or hand strength.



### MODEL 630-HS

Deluxe design: elastic side panels provide for superior cosmesis and comfort, operate independently of side pull straps. Zipper closure allows for easy application and removal. Casings for additional steels if desired.



### MODEL 622-HS

Similar in design to 630-HS, except with a 2" shorter back. Casings for additional steels if desired.

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	FRONT ADJUSTMENT	ELASTIC CLOSURE	HIP RELEASES	STEELS	EVEN DEV.	HIP SIZES
615-HS	Average to stocky	Average	FRONT 9.5" BACK 15.5"	White brocade	Three pull side lace	Hook-and-eye	At top and groin	1 pr. 13"	6"	30-48
615-HS	Average to stocky	Average	FRONT 9.5" BACK 15.5"	White brocade	Three pull side lace	Contact closure tabs	At top and groin	1 pr. 13"	6"	30-48
702-HS	Average to stocky	Average	FRONT 10" BACK 16"	White brocade	Three pull side lace	Hook-and-eye	Top, groin and bottom	1 pr. 13"	7"	30-48
630-HS	Average to full	Average to tall	FRONT 10" BACK 16"	Mesh/elastic	Three pull side lace	Zipper	Side panels all elastic	1 pr. 13"	9"	32-48
639-HS	Average to full	Average to tall	FRONT 10" BACK 16"	Mesh/elastic	Three pull side lace	Contact closure tabs	Side panels all elastic	1 pr. 13"	9"	32-48
622-HS	Average to full	Short to average	FRONT 10" BACK 14"	Mesh/elastic	Three pull side lace	Zipper	Side panels all elastic	1 pr. 11"	9"	32-48

## WOMEN'S LUMBOSACRAL SUPPORTS, LACED



### MODEL 674-HS

Generously cupped front for more effective control of a pendulous abdomen. Designed specifically for the heavy figure of average height.



### MODEL 664-HS

Hoke type garment is effective for control over muscular laxity in the torso - its broad size range allows for application even in adolescents.



### MODEL 370 SACROILIAC SUPPORT

Suitable for virtually all figure types. Includes removable sacral pad, adjustable lace back and single encircling pull strap. Hook-and-eye fastening in front. Elastic releases at groin.

Front depth 8", back depth 9"  
5" hip development  
Even hip sizes 28-44



### MODEL 676-HS

Similar in design to 674-HS, but more appropriate for the short, heavy figure.

MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	FRONT ADJUSTMENT	ELASTIC CLOSURE	HIP RELEASES	STEELS	EVEN DEV.	HIP SIZES
674-HS	Pendulous	Average	FRONT 13" BACK 17"	White brocade	Three pull side lace	Hook-and-eye	Top, groin, sides, back	1 pr. 13"	6"	36-56
676-HS	Pendulous	Average	FRONT 11.5" BACK 15"	White brocade	Three pull side lace	Hook-and-eye	Top, groin, sides, back	1 pr. 11"	7"	36-56
664-HS	Average	Average	FRONT 10.5" BACK 15.5"	White brocade	Pull strap/truss buckle	Hoke type	None	1 pr. 12"	9"	30-48

## WOMEN'S DORSOLUMBAR SUPPORTS, LACED



### MODEL 608-HS

Basic three-pull dorsolumbar with broad applications for the average to tall figure.



### MODEL 604-HS

Very similar in design to 608-HS, but with slightly shorter front and back depths.



### MODEL 606-HS

Similar to 604-HS, but with a shorter front. Cinch type shoulder straps with side front truss buckle fastening permits easier application by wearer with limited dexterity or hand strength.



### MODEL 605-HS

Another modification of 604-HS; cut from the same pattern, but with only a 6" hip development, and contact closure straps for the wearer with limited dexterity or hand strength.

### MODEL 775-HS

Specifically designed for the short, pendulous figure. Elastic side panels offer support and superior cosmesis for the pendulous abdomen. Side pulls operate independently of the elastic side panels.



MODEL	FIGURE	HEIGHT	DEPTH	MATERIAL	FRONT ADJUSTMENT	ELASTIC CLOSURE	HIP RELEASES	STEELS	EVEN DEV.	HIP SIZES
608-HS	Average	Average to tall	FRONT 10" BACK 24"	White brocade	Three pull side lace	Hook-and-eye	At top and groin	1 pr. 20"	8"	30-48
604-HS	Average	Average	FRONT 9" BACK 21"	White brocade	Three pull side lace	Hook-and-eye	At top and groin	1 pr. 18"	8"	30-48
606-HS	Average	Short to average	FRONT 8.5" BACK 21"	White brocade	Three pull side lace	Hook-and-eye	At top and groin	1 pr. 18"	8"	30-48
605-HS	Average/straight	Average	FRONT 9" BACK 21"	White brocade	Three pull side lace	Contact closure tabs	At top and groin	1 pr. 18"	6"	30-46
775-HS	Pendulous	Short to average	FRONT 8" BACK 18.5"	Brocade/elastic	Three pull side lace	Zipper	Side panels all elastic	1 pr. 14"	4"	30-48

## WOMEN'S ABDOMINAL AND MATERNITY SUPPORTS



### MODEL 591 TAYLOR (TLSO) INSERT

Converts a lumbosacral corset to a dorsolumbar by inserting the aluminum uprights into the stay casings and affixing the extension to the top of the corset (all attachments included). Removable and washable padded cover, foam shoulder straps.

Lengths: 17", 19", 21"



### MODEL 679 OBESITY SUPPORT

To support the pendulous abdomen; abdominal uplift is effected by two wide uplift panels. White brocade material with heavy elastic releases across the top and the groin. Left front hook-and-eye fastening. 13.5" reinforced cupped front; 17" lace back.

Even hip sizes 36-52



### MODEL 374 POST-OP SUPPORT

For use after abdominal surgery, or when abdominal muscles have lost their tone. Power net elastic sides and front, white brocade back. Back lace arrangement, bilateral pull straps. 9"-11" front, 7.5"-9.5" back  
Even hip sizes 30-48



### MODEL 927 RIB FRACTURE SUPPORT

Can be used as either a rib fracture support or post-operative binder. Elastic belt, 10" in back and 6" in front, with three vertical plastic stays (exterior flexible boning). Contact closure panel in front.

Sizes - Small: 29"-32" chest  
Medium: 33"-36" chest  
Large: 37"-40" chest  
X-Large: 41"-44" chest  
2X-Large: 45"-48" chest



### MODEL 683-HS MATERNITY SUPPORT

Helps reduce fatigue and backache which can accompany pregnancy. 11.5" woven elastic front with hook-and-eye opening exerts no pressure over abdomen; 16" back accommodates two 12" shapeable stays. Two pull side lacing adjustment; lacing is expandable; bottom should be pulled tighter than the top.

Even waist sizes 26-48



### MODEL 785-X MATERNITY SUPPORT

Two piece wraparound support with posterior pocket to accommodate thermoplastic support pad. Anterior panel exerts no pressure over the abdomen. Contact closure throughout. Application instructions included.

Sizes - S: 3-8 normal dress size  
M: 9-14 normal dress size  
L: 15-18 normal dress size  
XL: 20+ normal dress size