

SPRING TRUSSES



MODEL 531

Steel spring, rustproof, covered with 1" synthetic leather sleeve; 1.75" elastic band with adjustable keyhole buckle – infinite adjustment for maximum comfort. No. 394 inguinal pad, universal adjustment fixture, covered with synthetic leather. Plated steel shank with 1" long slot to allow for pad adjustment.

Even hip sizes 30-44, right and left



MODEL 348

Steel spring, rustproof, covered with 1" synthetic leather sleeve; 1.75" elastic band with adjustable keyhole buckle – infinite adjustment for maximum comfort. No. 391 scrotal pad, with perineal strap and universal adjustment fixture, covered with synthetic leather. Plated steel shank with 1" long slot to allow for pad adjustment.

Even hip sizes 30-44, right and left



MODEL 351

Steel spring, rustproof, covered with 1" synthetic leather sleeve; .875" leather fastening strap with stud holes spaced at 1" intervals. No. 394 inguinal pad, universal adjustment fixture, covered with synthetic leather. Plated steel shank with 1" long slot to allow for pad adjustment.

Even hip sizes 30-44, right and left



MODEL 301

Same as 531 above, but with the addition of 1" x 5.75" synthetic leather back pad (shown above) to minimize pressure over lower back musculature.

Even hip sizes 30-44, right and left



MODEL 348-K

Same as 348 above, but with the addition of 1" x 5.75" synthetic leather back pad (shown above) to minimize pressure over lower back musculature.

Even hip sizes 30-44, right and left



MODEL 302

Steel spring, rustproof, covered with 1" synthetic leather sleeve; .875" leather fastening strap with stud holes spaced at 1" intervals. No. 391 scrotal pad, with perineal strap and universal adjustment fixture, covered with synthetic leather. Plated steel shank with 1" long slot to allow for pad adjustment.

Even hip sizes 30-44, right and left

HOOD AND DOUBLE SPRING, ELASTIC TRUSSES



MODEL 29-T

Elastic band, 1.75" wide, with adjustable keyhole buckle – infinite adjustment for maximum comfort. Scrotal pad – sewn onto the elastic band – with elastic perineal strap featuring claw buckle posterior fastening. Pad is nonadjustable, covered with synthetic leather.

Even hip sizes 30-44, right and left



MODEL 352

Two (right and left) rustproof steel springs covered with 1" synthetic leather sleeves, joined together by no. 198 tapered back pad, 7" x 2.5" at apex, which allows for circumferential adjustment and lateral pad placement. Two no 394 inguinal pads, with universal adjustment fixture and synthetic leather covering. Two (no. 59) plated steel shanks. No. 58 .75" anterior leather fastening strap, 10" long, with stud holes spaced at 1" intervals.

Even hip sizes 30-44

MODEL 344

Two (right and left) rustproof steel springs covered with 1" synthetic leather sleeves, flared at distal ends to minimize discomfort over lower back musculature; plated steel anterior connector bar with 1.5" long slots for lateral pad placement. Two no 389 inguinal pads, with universal adjustment fixture and synthetic leather covering. No. 58 .75" posterior leather fastening strap, 10" long, with stud holes spaced at 1" intervals.

Even hip sizes 30-44



MODEL 30-T

Elastic band, 1.75" wide, with adjustable keyhole buckle – infinite adjustment for maximum comfort. Right and left scrotal pads, universal adjustment fixture, with elastic perineal straps, featuring claw buckle posterior fastening. Plated anterior bar, 6" long, with 1" slots for lateral pad adjustment.

Even hip sizes 30-44



MODEL 331

Same as 352 above, except that two no. 391 scrotal pads, with perineal strap and synthetic leather covering, are utilized.

Even hip sizes 30-44



MODEL 247-F

1.75" wide elastic band – fastens to anterior flange with adjustable keyhole buckle, and over pad with double .875" web loop arrangement and snap billet hooks. Oval inguinal pad, covered with synthetic leather, adjusts to any angle. Adjustable elastic perineal strap with claw buckle posterior fastening – attaches to anterior flange beneath lower snap billet hook.

Even hip sizes 30-44, right and left



MODEL 248-F

As 247-F above, but with both right and left inguinal pads and differently shaped anterior flange. For use with bilateral hernias.

Even hip sizes 30-44

TRUSS PARTS AND PADS



MODEL 312

Same as 344 on page 32, except that one inguinal pad is replaced with a false pad (shown on right) in cases of single hernia for placement over the unaffected side.

Even hip sizes 30-44



MODEL 198 BACK PAD

For double spring truss; covered with synthetic leather, with holes for lateral spring adjustment; 7" x 2.5" at apex.

MODEL 197 BACK PAD

For single spring truss; covered with synthetic leather; 5.75" x 2.5".

MODEL 58 LEATHER CONNECTING STRAP

For hood or double spring truss; 10" x .75" with stud holes spaced at 1" intervals.

MODEL 67 FRONT BAR FOR HOOD TRUSS

Chrome plated, 9.5" long with 1.5" long slots for pad adjustment.

MODEL 55 ELASTIC PERINEAL STRAP

For use on Boston models; 20" long, with snap billet hook fastener.

MODEL 56 NON-ELASTIC PERINEAL STRAP

For use on Boston models; 20" long, moleskin and non-elastic web.

MODEL 59 SHANK

Chrome plated, with 1" long slot for pad adjustment; specify right or left.

MODEL 556 PERINEAL STRAP FASTENER

For elastic trusses; non-elastic web loop with claw buckle fastener.

MODEL 445 elastic web - 1.75" wide

MODEL 445-C elastic web - .875" wide

MODEL 445-G non-elastic web - .875" wide

MODEL 445-H non-elastic web - .75" wide

MODEL 516 LEATHER STRAP FASTENER

For spring trusses; leather loop with claw buckle fastener.

MODEL 389 INGUINAL PAD

For use on hood trusses; oval shape, covered with synthetic leather; 3" x 2.25" x 1.25" thick.

MODEL 388 INGUINAL PAD

Oval shape, covered with synthetic leather; universal fixture. 3" x 2.25" x 1.25" thick.

MODEL 394 FRENCH SHAPE PAD

Full bottom for comfort. Covered with synthetic leather; universal fixture. 3.5" x 2.25" x 1" thick.

MODEL 391 SCROTAL PAD

Includes moleskin/leather perineal strap and no. 516 fastener. Covered with synthetic leather. 5" x 3.25" x 1.25" thick; specify right or left.

TRUSS PARTS AND PADS



MODEL 35-P UMBILICAL TRUSS

Double 1.75" elastic band – 3.5" total width – with double non-elastic web pull strap adjustment and truss buckle fastening. Spnge rubber pad – 4" x 6" – with center bump.

Even hip sizes 30-44



MODEL 174 UMBILICAL TRUSS

2" white elastic web with hook-and-loop fastening. Sponge rubber pad – 2.5" x 2.5" – with center bump.

Even hip sizes 10-24



MODEL 490-4E ATHLETIC SUPPORT

Stretchy gray knitted elastic, 9" deep all around. Flexible spiral boning in front and back, minimizes rolling and migration. Detachable knitted gray fly pouch with leg straps – snap fastening (no. 490-4PE).

Sizes – X-Small: 27"-30" waist
 Small: 31"-34" waist
 Medium: 35"-38" waist
 Med. Large: 39"-42" waist
 Large: 43"-46" waist
 X-Large: 47"-50" waist



MODEL 243 SUSPENSORY

1.75" elastic web waistband with claw buckle fastening. Flannel lined mesh pouch – .75" elastic perineal straps with claw buckle fastening. Sized according to pouch. Medium is the size most often selected.

Sizes – Medium
 Large
 X-Large
 2X-Large



MODEL 270 SUSPENSORY

Soft mesh pouch, spaced to provide comfortable support with the use of under- straps. Non-elastic web waistband, attached to ensure correct positioning of pouch. Sized according to pouch. Medium is the size most often selected.

Sizes – X-Small
 Small
 Medium
 Large
 X-Large
 2X-Large
 3X-Large



MODEL 272 SUSPENSORY

Same as model 270 at left, but with the addition of leg straps for additional control of pouch position.

Sizes – X-Small
 Small
 Medium
 Large
 X-Large
 2X-Large
 3X-Large